

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP08724 - IMMUNIZING PEPTIDE

SKU: EBP08724

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μα

Shipped at ambient temperature, store at -20°C.

NP_004296.1; NP_647593.1; NP_647594.1; NP_647595.1; NP_647596.1; NP_647597.1; NP_647598.1; Accession

NP_647599.1; NP_647600.1; NP_647601.1 ID

Peptide with sequence C-SDNAPAKGNKSPS, from the internal region of the protein sequence according to

Immunogen NP 004296.1; NP 647593.1; NP 647594.1; NP 647595.1; NP 647596.1; NP 647597.1; NP 647598.1;

NP 647599.1; NP 647600.1; NP 647601.1.

The peptide was used in the production of EB08724 - Goat Anti-BIN1 Antibody. It is currently untested in **Product**

Comments blocking applications.

Peptide

C-SDNAPAKGNKSPS Sequence

Shipping

Instructions Ambient

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade