

# EBP08734 - IMMUNIZING PEPTIDE

**SKU:** EBP08734



---

## SPECIFICATIONS

|                              |   |
|------------------------------|---|
| <b>Formulation</b>           | 100µg of Lyophilised peptide, formulation De-ionized (DI) water.  |
| <b>Unit Size</b>             | 100 µg  |
| <b>Storage Instructions</b>  | Shipped at ambient temperature, store at -20°C.   |
| <b>Accession ID</b>          | NP_004987.2; NP_063915.2; NP_063953.2; NP_063954.2; NP_063955.2   |
| <b>Immunogen</b>             | Peptide with sequence C-HQSDLKVDENQKAYY, from the internal region of the protein sequence according to NP_004987.2; NP_063915.2; NP_063953.2; NP_063954.2; NP_063955.2.   |
| <b>Product Comments</b>      | The peptide was used in the production of EB08734 - Goat Anti-ABCC1 Antibody. It is currently untested in blocking applications.  |
| <b>Peptide Sequence</b>      | C-HQSDLKVDENQKAYY   |
| <b>Shipping Instructions</b> | Ambient   |
| <b>Product Grade</b>         | <a href="https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png">https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png</a> |

## DOCUMENTS

- [Data Sheet](#)