

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP08740 - IMMUNIZING PEPTIDE**

**SKU:** EBP08740

## **SPECIFICATIONS**

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

**Unit Size** 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

**Accession** 

 $NP\_058198.2;\ NP\_997057.1;\ NP\_004906.3;\ NP\_997510.1;\ NP\_997511.1;\ NP\_997512.1$ ID

 $\textbf{Immunogen} \begin{array}{l} \textbf{Peptide with sequence C-YFVLRYKIRAER, from the C Terminus of the protein sequence according to NP_058198.2;} \\ \textbf{NP_997057.1; NP_004906.3; NP_997510.1; NP_997511.1; NP_997512.1.} \end{array}$ 

The peptide was used in the production of EB08740 - Goat Anti-ABCG1 Antibody. It is currently untested in blocking **Product** 

Comments applications.

**Peptide** 

C-YFVLRYKIRAER Sequence

Shipping

Ambient Instructions

**Product** 

 $https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_plus\_medium.png$ Grade