



EBP08740 - IMMUNIZING PEPTIDE

SKU: EBP08740



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_058198.2; NP_997057.1; NP_004906.3; NP_997510.1; NP_997511.1; NP_997512.1
Immunogen	Peptide with sequence C-YFVLRYSKIRAER, from the C Terminus of the protein sequence according to NP_058198.2; NP_997057.1; NP_004906.3; NP_997510.1; NP_997511.1; NP_997512.1.
Product Comments	The peptide was used in the production of EBP08740 - Goat Anti-ABCG1 Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-YFVLRYSKIRAER
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png