

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP08777 - IMMUNIZING PEPTIDE

SKU: EBP08777

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage Shipped at ambient temperature, store at -20°C.

ID

Accession NP_000925.2; NP_001159393.1

The peptide was used in the production of EB08777 - Goat Anti-Alpha2-antiplasmin Antibody. It is currently untested **Product**

Comments in blocking applications.

Peptide

C-KDFLQSLKGFPRGDK Sequence

Shipping

Ambient Instructions

Product

 $https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png$ Grade