



# EBP08796 - IMMUNIZING PEPTIDE

**SKU:** EBP08796



---

## SPECIFICATIONS

**Formulation** 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

**Unit Size** 100 µg

**Storage Instructions** Shipped at ambient temperature, store at -20°C.

**Accession ID** NP\_057263.1; NP\_003297.1; NP\_001129427.1; NP\_001129429.1; NP\_001129428.1; NP\_001129430.1

**Immunogen** Peptide with sequence C-KHAKEEDSSIDYD, from the internal region (near C Terminus) of the protein sequence according to NP\_057263.1; NP\_003297.1; NP\_001129427.1; NP\_001129429.1; NP\_001129428.1; NP\_001129430.1.

**Product Comments** The peptide was used in the production of EB08796 - Goat Anti-TRPC4 Antibody. It is currently untested in blocking applications.

**Peptide Sequence** C-KHAKEEDSSIDYD

**Shipping Instructions** Ambient

**Product Grade** [https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\\_medium.png](https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png)

## DOCUMENTS

- [Data Sheet](#)