

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP08796 - IMMUNIZING PEPTIDE

SKU: EBP08796

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μα

Shipped at ambient temperature, store at -20°C.

Accession

NP 057263.1; NP 003297.1; NP 001129427.1; NP 001129429.1; NP 001129428.1; NP 001129430.1 ID

Peptide with sequence C-KHAKEEDSSIDYD, from the internal region (near C Terminus) of the protein sequence

Immunogen according to NP 057263.1; NP 003297.1; NP 001129427.1; NP 001129429.1; NP 001129428.1;

NP 001129430.1.

The peptide was used in the production of EB08796 - Goat Anti-TRPC4 Antibody. It is currently untested in **Product**

Comments blocking applications.

Peptide

C-KHAKEEDSSIDYD Sequence

Shipping

Instructions Ambient

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade