

EBP08830 - IMMUNIZING PEPTIDE

SKU: EBP08830



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	
Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_002551.2; NP_001243725.1; NP_001248326.1
Immunogen	Peptide with sequence C-YREKKYKYVEDYEQ, from the C Terminus of the protein sequence according to NP_002551.2; NP_001243725.1; NP_001248326.1.
Product	The peptide was used in the production of EBP08830 - Goat Anti-P2X4 / P2X4R Antibody. It is currently untested
Comments	in blocking applications.
Peptide Sequence	C-YREKKYKYVEDYEQ
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)