

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP08834 - IMMUNIZING PEPTIDE

SKU: EBP08834

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

100 μg **Unit Size**

Storage

Shipped at ambient temperature, store at -20°C. Instructions

Accession

NP 443134.2 ID

Peptide with sequence C-SRGTPNRERKQGEQS, from the internal region of the protein sequence according to **Immunogen**

NP 443134.2.

The peptide was used in the production of EB08834 - Goat Anti-LIP1 / LKB1IP Antibody. It is currently untested **Product**

Comments in blocking applications.

Peptide

C-SRGTPNRERKQGEQS Sequence

Shipping

Instructions Ambient

Product

 $https://prod-vector-labs-pimcore-assets.s 3. us-east-1. amazonaws.com/assets/products/image/elite_medium.png$ Grade