

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP08847 - IMMUNIZING PEPTIDE**

**SKU:** EBP08847

## **SPECIFICATIONS**

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$ 

100 μg **Unit Size** 

Storage

Shipped at ambient temperature, store at -20°C. Instructions

**Accession** 

ID

NP 071407.3; NP 443068.1

Peptide with sequence C-YNISLYENVTVGTS , from the internal region of the protein sequence according to Immunogen

NP 071407.3; NP 443068.1.

The peptide was used in the production of EB08847 - Goat Anti-CDH23 / USH1D Antibody. It is currently **Product** 

**Comments** untested in blocking applications.

**Peptide** 

C-YNISLYENVTVGTS Sequence

Shipping

Instructions Ambient

**Product** Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png