

EBP08878 - IMMUNIZING PEPTIDE

SKU: EBP08878



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_001824.1; NP_009027.1; NP_001070145.1; NP_001825.1; NP_009028.1
Immunogen	Peptide with sequence CDFNPKSSQAKD, from the internal region of the protein sequence according to NP_001824.1; NP_009027.1; NP_001070145.1; NP_001825.1; NP_009028.1.
Product Comments	The peptide was used in the production of EBP08878 - Goat Anti-CTLA + CTLB Antibody. It is currently untested in blocking applications.
Peptide Sequence	CDFNPKSSQAKD
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)