

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP08890 - IMMUNIZING PEPTIDE

SKU: EBP08890

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_005328.2; NP_001171905.1 ID

Immunogen Peptide with sequence C-RHTENVTKTKTPE, from the internal region of the protein sequence according to NP_005328.2; NP_001171905.1.

The peptide was used in the production of EB08890 - Goat Anti-HEM1 Antibody. It is currently untested in blocking **Product**

Comments applications.

Peptide

C-RHTENVTKTKTPE

Sequence Shipping

Ambient Instructions

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade