

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP08892 - IMMUNIZING PEPTIDE

SKU: EBP08892

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

ID

NP_758438.2

Immunogen Peptide with sequence C-HHQSFHYEFQHSHP, from the internal region of the protein sequence according to NP_758438.2.

Product The peptide was used in the production of EB08892 - Goat Anti-TFAP2D Antibody. It is currently untested in blocking

Comments applications.

Peptide

Sequence

C-HHQSFHYEFQHSHP

Shipping

Ambient Instructions

Product

Grade

 $https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png$