

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP08897 - IMMUNIZING PEPTIDE

SKU: EBP08897

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μα

Shipped at ambient temperature, store at -20°C.

Accession NP_002561.1; NP_612192.1; NP_612195.1; NP_612197.1; NP_612196.1; NP_612198.2; NP_612193.1;

ID NP_612194.1

Peptide with sequence C-YHTFSAHQSRSR, from the internal region of the protein sequence according to

Immunogen NP 002561.1; NP 612192.1; NP 612195.1; NP 612197.1; NP 612196.1; NP 612198.2; NP 612193.1;

The peptide was used in the production of EB08897 - Goat Anti-PCSK6 Antibody. It is currently untested in **Product**

Comments blocking applications.

Peptide

C-YHTFSAHQSRSR Sequence

Shipping

Instructions Ambient

Product https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade