

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP08931 - IMMUNIZING PEPTIDE

SKU: EBP08931

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Shipped at ambient temperature, store at -20°C.

Accession

 $NP_001008491.1; \ NP_001269901.1; \ NP_001269902.1; \ NP_001307958.1; \ NP_001336231.1$ ID

Peptide with sequence SKQQPTQFINPET-C, from the N terminus of the protein sequence according to Immunogen NP_001008491.1; NP_001269901.1; NP_001269902.1; NP_001307958.1; NP_001336231.1.

The peptide was used in the production of EB08931 - Goat Anti-SEPT2 Antibody. It is currently untested in blocking **Product**

Comments applications.

Peptide

SKQQPTQFINPET-C Sequence

Shipping

Ambient Instructions

Product

 $https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png$ Grade