

EBP08932-B - IMMUNIZING PEPTIDE

SKU: EBP08932-B



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	Shipped at ambient temperature, store at -20°C.
Instructions	
Accession ID	NP_665798.1; NP_055944.2; NP_665801.1
Immunogen	Peptide with sequence C-DEVNAFKQRKTA., from the internal region of the protein sequence according to NP_665798.1; NP_055944.2; NP_665801.1.
Product Comments	The peptide was used in the production of EBP08932-B - Goat Anti-SEPT6, Biotinylated Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-DEVNAFKQRKTA.
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)