

EBP08982 - IMMUNIZING PEPTIDE

SKU: EBP08982



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	
Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_060405.3; NP_001077376.1; NP_001120979.1;
Immunogen	Peptide with sequence C-ESRKTYGQSADGAE, from the internal region of the protein sequence according to NP_060405.3; NP_001077376.1; NP_001120979.1;.
Product	The peptide was used in the production of EB08982 - Goat Anti-BANK1 Antibody. It is currently untested in
Comments	blocking applications.
Peptide Sequence	C-ESRKTYGQSADGAE
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png