

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP09008 - IMMUNIZING PEPTIDE

SKU: EBP09008

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μα

Storage

Shipped at ambient temperature, store at -20°C. Instructions

Accession

NP_067023.1; NP_001123506.1; NP_001123505.1; NP_001123507.1; NP_001123508.1; ID

Peptide with sequence PPDSALESPFEE-C, from the N Terminus of the protein sequence according to Immunogen

NP 067023.1; NP 001123506.1; NP 001123505.1; NP 001123507.1; NP 001123508.1;.

The peptide was used in the production of EB09008 - Goat Anti-PHR1 Antibody. It is currently untested in **Product**

Comments blocking applications. **Peptide**

Sequence

PPDSALESPFEE-C

Shipping

Instructions Ambient

Product Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png