

EBP09016 - IMMUNIZING PEPTIDE

SKU: EBP09016



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	Shipped at ambient temperature, store at -20°C.
Instructions	
Accession ID	NP_005092.1
Immunogen	Peptide with sequence C-RNNKGRSSTYEVR, from the internal region of the protein sequence according to NP_005092.1.
Product	The peptide was used in the production of EB09016 - Goat Anti-ISG15 Antibody. It is currently untested in blocking applications.
Comments	
Peptide Sequence	C-RNNKGRSSTYEVR
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png