

EBP09093 - IMMUNIZING PEPTIDE

SKU: EBP09093



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	
Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_001461.1; NP_068703.1; NP_068704.2; NP_068705.2
Immunogen	Peptide with sequence DRLSCDGSRVHLLYK, from the C Terminus of the protein sequence according to NP_001461.1; NP_068703.1; NP_068704.2; NP_068705.2.
Product	The peptide was used in the production of EB09093 - Goat Anti-GABBR1 Antibody. It is currently untested in
Comments	blocking applications.
Peptide Sequence	DRLSCDGSRVHLLYK
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)