

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP09132 - IMMUNIZING PEPTIDE

SKU: EBP09132

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_945350.1; NP_996766.1; NP_996767.1; NP_996768.1 ID

Immunogen Peptide with sequence ERLHTKDKVPDSIAD -C, from the N Terminus (near) of the protein sequence according to NP_945350.1; NP_996766.1; NP_996767.1; NP_996768.1.

The peptide was used in the production of EB09132 - Goat Anti-SLC26A5 / prestin Antibody. It is currently untested Product

Comments in blocking applications.

Peptide

ERLHTKDKVPDSIAD-C Sequence

Instructions Ambient Shipping

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade