

EBP09132 - IMMUNIZING PEPTIDE

SKU: EBP09132



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_945350.1; NP_996766.1; NP_996767.1; NP_996768.1
Immunogen	Peptide with sequence ERLHTKDKVPDSIAD -C, from the N Terminus (near) of the protein sequence according to NP_945350.1; NP_996766.1; NP_996767.1; NP_996768.1.
Product Comments	The peptide was used in the production of EBP09132 - Goat Anti-SLC26A5 / prestin Antibody. It is currently untested in blocking applications.
Peptide Sequence	ERLHTKDKVPDSIAD -C
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

DOCUMENTS

- [Data Sheet](#)