

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP09140 - IMMUNIZING PEPTIDE

SKU: EBP09140

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_067042.1 ID

 $\begin{tabular}{ll} \textbf{Immunogen} & Peptide with sequence C-NPRLGSQSTNSSYT , from the internal region of the protein sequence according to NP_067042.1. \end{tabular}$

The peptide was used in the production of EB09140 - Goat Anti-JAM2 / JAM-B Antibody. It is currently untested in **Product**

Comments blocking applications.

Peptide

C-NPRLGSQSTNSSYT Sequence

Shipping

Ambient Instructions

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade