

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP09152 - IMMUNIZING PEPTIDE

SKU: EBP09152

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Shipped at ambient temperature, store at -20°C.

Accession

NP_003237.2 ID

 $\label{eq:local_problem} \textbf{Immunogen} \begin{tabular}{l} \textbf{Peptide with sequence NRIPESGGDNSVFD-C, from the N Terminus of the protein sequence according to NP_003237.2.} \end{tabular}$

Product The peptide was used in the production of EB09152 - Goat Anti-Thrombospondin Antibody. It is currently untested in

Comments blocking applications.

Peptide

NRIPESGGDNSVFD-C Sequence

Shipping

Ambient Instructions

Product

 $https://prod-vector-labs-pimcore-assets.s 3. us-east-1. amazonaws. com/assets/products/image/elite_plus_medium.png$ Grade