

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP09167-B - IMMUNIZING PEPTIDE

SKU: EBP09167-B

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μα

Storage

Shipped at ambient temperature, store at -20°C. Instructions

Accession

NP 053733.2; NP 005234.1; NP 001156757.1; NP 001156758.1; NP 001156759.1 ID

Peptide with sequence C-TSYDQSSYSQQNTYG., from the internal region of the protein sequence according to Immunogen

NP 053733.2; NP 005234.1; NP 001156757.1; NP 001156758.1; NP 001156759.1.

The peptide was used in the production of EB09167-B - Goat Anti-EWS / EWSR1, Biotinylated Antibody. It is **Product**

Comments currently untested in blocking applications.

Peptide

C-TSYDQSSYSQQNTYG. Sequence

Shipping

Instructions Ambient

Product Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png