

EBP09167-B - IMMUNIZING PEPTIDE

SKU: EBP09167-B



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_053733.2; NP_005234.1; NP_001156757.1; NP_001156758.1; NP_001156759.1
Immunogen	Peptide with sequence C-TSYDQSSYSQQNTYG., from the internal region of the protein sequence according to NP_053733.2; NP_005234.1; NP_001156757.1; NP_001156758.1; NP_001156759.1.
Product Comments	The peptide was used in the production of EB09167-B - Goat Anti-EWS / EWSR1, Biotinylated Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-TSYDQSSYSQQNTYG.
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)