

EBP09167 - IMMUNIZING PEPTIDE

SKU: EBP09167



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_053733.2; NP_005234.1; NP_001156757.1; NP_001156758.1; NP_001156759.1
Immunogen	Peptide with sequence C-TSYDQSSYSQQNTYG, from the internal region of the protein sequence according to NP_053733.2; NP_005234.1; NP_001156757.1; NP_001156758.1; NP_001156759.1.
Product Comments	The peptide was used in the production of EB09167 - Goat Anti-EWS / EWSR1 Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-TSYDQSSYSQQNTYG
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)