

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP09170 - IMMUNIZING PEPTIDE**

**SKU:** EBP09170

## **SPECIFICATIONS**

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$ 

**Unit Size** 100 μα

Storage

Shipped at ambient temperature, store at -20°C. Instructions

**Accession** 

NP 490597.1; NP 055114.1 ID

Peptide with sequence C-PETSNPNKPKRQTNQ, from the internal region (near N Terminus) of the protein **Immunogen** 

sequence according to NP 490597.1; NP 055114.1.

The peptide was used in the production of EB09170 - Goat Anti-BRD4 Antibody. It is currently untested in **Product** 

**Comments** blocking applications. **Peptide** C-PETSNPNKPKRQTNQ

Sequence

Shipping Instructions Ambient

**Product** https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png Grade