



EBP09172 - IMMUNIZING PEPTIDE

SKU: EBP09172



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	Shipped at ambient temperature, store at -20°C.
Instructions	
Accession ID	NP_003880.3; NP_071285.1; NP_148934.1
Immunogen	Peptide with sequence EQFAITLKSIECNR, from the internal region of the protein sequence according to NP_003880.3; NP_071285.1; NP_148934.1.
Product Comments	The peptide was used in the production of EB09172 - Goat Anti-PXR / NR1I2 Antibody. It is currently untested in blocking applications.
Peptide Sequence	EQFAITLKSIECNR
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)