

EBP09173 - IMMUNIZING PEPTIDE

SKU: EBP09173



SPECIFICATIONS

| | |
|------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Formulation | 100µg of Lyophilised peptide, formulation De-ionized (DI) water. |
| Unit Size | 100 µg |
| Storage | Shipped at ambient temperature, store at -20°C. |
| Instructions | |
| Accession ID | NP_001003891.1; NP_056973.2 |
| Immunogen | Peptide with sequence C-PPQSQTQPLVSQAQ, from the internal region of the protein sequence according to NP_001003891.1; NP_056973.2. |
| Product Comments | The peptide was used in the production of EB09173 - Goat Anti-PCQAP / MED15 Antibody. It is currently untested in blocking applications. |
| Peptide Sequence | C-PPQSQTQPLVSQAQ |
| Shipping Instructions | Ambient |
| Product Grade | https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png |

DOCUMENTS

- [Data Sheet](#)