

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP09186 - IMMUNIZING PEPTIDE

SKU: EBP09186

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_003950.1 ID

Immunogen Peptide with sequence C-RQDHQLDKKDGIYP, from the C Terminus of the protein sequence according to NP_003950.1.

The peptide was used in the production of EB09186 - Goat Anti-HIP1R/ HIP12 Antibody. It is currently untested in **Product**

Comments blocking applications. **Peptide**

C-RQDHQLDKKDGIYP Sequence

Shipping

Ambient Instructions

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade