

EBP09195 - IMMUNIZING PEPTIDE

SKU: EBP09195



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_002986.1; NP_003166.1
Immunogen	Peptide with sequence C-RDVVRSMRKSNTRN, from the internal region of the protein sequence according to NP_002986.1; NP_003166.1.
Product Comments	The peptide was used in the production of EBP09195 - Goat Anti-XCL1 / XCL2 Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-RDVVRSMRKSNTRN
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspising_medium.png

DOCUMENTS

- [Data Sheet](#)