

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP09195 - IMMUNIZING PEPTIDE**

**SKU:** EBP09195

## **SPECIFICATIONS**

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

**Unit Size** 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP\_002986.1; NP\_003166.1 ID

 $\label{eq:lower_lower} \textbf{Immunogen} \ \, \frac{\text{Peptide with sequence C-RDVVRSMDRKSNTRN, from the internal region of the protein sequence according to NP_002986.1; NP_003166.1.}$ 

The peptide was used in the production of EB09195 - Goat Anti-XCL1 / XCL2 Antibody. It is currently untested in **Product** 

**Comments** blocking applications.

Peptide

C-RDVVRSMDRKSNTRN

Sequence

Ambient Instructions Shipping

**Product** 

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring\_medium.png Grade