

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP09196 - IMMUNIZING PEPTIDE

SKU: EBP09196

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_005742.4; NP_671714.1 ID

 $\label{eq:lower_lower_lower} \textbf{Immunogen} \ \ \frac{\text{Peptide with sequence QRKAQSDGQSPSKK-C, from the N Terminus of the protein sequence according to NP_005742.4; NP_671714.1.}$

The peptide was used in the production of EB09196 - Goat Anti-AKAP9 / YOTIAO Antibody. It is currently untested Product

Comments in blocking applications.

Peptide

QRKAQSDGQSPSKK-C Sequence

Shipping

Ambient Instructions

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade