

EBP09202 - IMMUNIZING PEPTIDE

SKU: EBP09202



SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 µg

Storage Instructions Shipped at ambient temperature, store at -20°C.

Accession ID NP_443095.2

Immunogen Peptide with sequence NHLIEGSQKCVAE, from the internal of the protein sequence according to NP_443095.2.

Product Comments The peptide was used in the production of EB09202 - Goat Anti-SCGB3A1 / UGRP2 Antibody. It is currently untested in blocking applications.

Peptide Sequence NHLIEGSQKCVAE

Shipping Instructions Ambient

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

DOCUMENTS

- [Data Sheet](#)