

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP09202 - IMMUNIZING PEPTIDE

SKU: EBP09202

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_443095.2 ID

Immunogen Peptide with sequence NHLIEGSQKCVAE, from the internal of the protein sequence according to NP_443095.2.

Product The peptide was used in the production of EB09202 - Goat Anti-SCGB3A1 / UGRP2 Antibody. It is currently

Comments untested in blocking applications.

Peptide

Sequence

NHLIEGSQKCVAE

Shipping Ambient Instructions

Product

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png