

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP09212 - IMMUNIZING PEPTIDE**

**SKU:** EBP09212

## **SPECIFICATIONS**

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

**Unit Size** 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

ID

NP\_062816.2; NP\_001153602.1; NP\_001153604.1; NP\_1153605.1; NP\_001153606.1

Peptide with sequence C-ESTDALSLPHVKLK, from the C Terminus of the protein sequence according to

NP\_062816.2; NP\_001153602.1; NP\_001153604.1; NP\_1153605.1; NP\_001153606.1.

The peptide was used in the production of EB09212 - Goat Anti-KCNQ5 Antibody. It is currently untested in Product

Comments blocking applications.

Peptide

C-ESTDALSLPHVKLK Sequence

Shipping

Ambient Instructions

**Product** 

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring\_medium.png Grade