

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP09217 - IMMUNIZING PEPTIDE**

**SKU:** EBP09217

## **SPECIFICATIONS**

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

**Unit Size** 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP\_734465.2 ID

Peptide with sequence C-EALQQDDGDLN, from the internal region (near C-Terminus) of the protein sequence

according to NP\_734465.2.

The peptide was used in the production of EB09217 - Goat Anti-UNC5B Antibody. It is currently untested in **Product** 

**Comments** blocking applications.

Peptide

C-EALQQDDGDLN Sequence

Ambient Instructions Shipping

**Product** 

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring\_medium.png Grade