

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP09223 - IMMUNIZING PEPTIDE

SKU: EBP09223

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_116274.3; NP_001268433.1 ID

Peptide with sequence CSQLAQPTLRLPRT, from the internal region (near C terminus) of the protein sequence

according to NP_116274.3; NP_001268433.1.

The peptide was used in the production of EB09223 - Goat Anti-ATG4D Antibody. It is currently untested in blocking **Product**

Comments applications.

Peptide

CSQLAQPTLRLPRT Sequence

Shipping

Ambient Instructions

Product

 $https://prod-vector-labs-pimcore-assets.s 3. us-east-1. amazonaws. com/assets/products/image/elite_plus_medium.png$ Grade