

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP09243 - IMMUNIZING PEPTIDE

SKU: EBP09243

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

100 μg **Unit Size**

Storage

Shipped at ambient temperature, store at -20°C. Instructions

Accession

NP 000835.1; NP 870989.1 ID

Peptide with sequence NCKLTISGSKKEDT , from the internal region of the protein sequence according to Immunogen

NP 000835.1; NP 870989.1.

The peptide was used in the production of EB09243 - Goat Anti-GRM7 Antibody. It is currently untested in **Product**

Comments blocking applications. **Peptide**

Sequence

NCKLTISGSKKEDT

Shipping

Instructions Ambient

Product Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png