

EBP09248 - IMMUNIZING PEPTIDE

SKU: EBP09248



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_005333.2
Immunogen	Peptide with sequence C-KFDGAKGPAKVARKK, from the internal region (near C Terminus) of the protein sequence according to NP_005333.2.
Product Comments	The peptide was used in the production of EB09248 - Goat Anti-HMGB3 / HMG4 Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-KFDGAKGPAKVARKK
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)