

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP09285 - IMMUNIZING PEPTIDE

SKU: EBP09285

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_005609.3 ID

 $\begin{tabular}{l|l} \textbf{Immunogen} \\ \hline & \textbf{Peptide with sequence C-ATQRHLTVGEEWSQD, from the internal region of the protein sequence according to NP_005609.3.} \\ \hline \end{tabular}$

Product The peptide was used in the production of EB09285 - Goat Anti-DLL1 Antibody. It is currently untested in blocking

Comments applications.

Peptide

C-ATQRHLTVGEEWSQD Sequence

Shipping

Ambient Instructions

Product

 $https://prod-vector-labs-pimcore-assets.s 3. us-east-1. amazonaws. com/assets/products/image/elite_plus_medium.png$ Grade