

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP09301 - IMMUNIZING PEPTIDE

SKU: EBP09301

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage Shipped at ambient temperature, store at -20°C.

Accession

NP_000156.1 ID

Product The peptide was used in the production of EB09301 - Goat Anti-Connexin 43 / GJA1 Antibody. It is currently untested

Comments in blocking applications.

Peptide

C-QPFDFPDDNQNSKK Sequence

Shipping

Ambient Instructions

Product

 $https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png$ Grade