

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP09345 - IMMUNIZING PEPTIDE

SKU: EBP09345

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_116190.2; NP_001192258.1 ID

Peptide with sequence C-KQDGESYKNPGKPD, from the internal region (near C-Terminus) of the protein sequence

according to NP_116190.2; NP_001192258.1.

The peptide was used in the production of EB09345 - Goat Anti-JAMC Antibody. It is currently untested in blocking Product

Comments applications.

Peptide

C-KQDGESYKNPGKPD Sequence

Shipping

Ambient Instructions

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade