

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP09379 - IMMUNIZING PEPTIDE

SKU: EBP09379

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μα

Storage

Shipped at ambient temperature, store at -20°C. Instructions

Accession

NP 003799.1; NP 742084.1; NP 742085.1 ID

Peptide with sequence C-PRARAKLNLSPHGT, from the C-Terminus (near) of the protein sequence according to Immunogen

NP 003799.1; NP 742084.1; NP 742085.1.

The peptide was used in the production of EB09379 - Goat Anti-TNFSF13 / APRIL Antibody. It is currently **Product**

Comments untested in blocking applications.

Peptide

C-PRARAKLNLSPHGT Sequence

Shipping

Instructions Ambient

Product Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png