

EBP09379 - IMMUNIZING PEPTIDE

SKU: EBP09379



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	Shipped at ambient temperature, store at -20°C.
Instructions	
Accession ID	NP_003799.1; NP_742084.1; NP_742085.1
Immunogen	Peptide with sequence C-PRARAKLNLSPHGT, from the C-Terminus (near) of the protein sequence according to NP_003799.1; NP_742084.1; NP_742085.1.
Product Comments	The peptide was used in the production of EB09379 - Goat Anti-TNFSF13 / APRIL Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-PRARAKLNLSPHGT
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)