

EBP09387 - IMMUNIZING PEPTIDE

SKU: EBP09387



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	
Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_874365.3; NP_056171.3
Immunogen	Peptide with sequence C-PEGPGKEKELPGQ, from the internal region of the protein sequence according to NP_874365.3; NP_056171.3.
Product	The peptide was used in the production of EB09387 - Goat Anti-SCRIB / Scribble Antibody. It is currently
Comments	untested in blocking applications.
Peptide Sequence	C-PEGPGKEKELPGQ
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png