

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP09403 - IMMUNIZING PEPTIDE**

**SKU:** EBP09403

## **SPECIFICATIONS**

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$ 

**Unit Size** 100 μα

Storage

Shipped at ambient temperature, store at -20°C. Instructions

**Accession** 

NP 005989.3; NP 001008883.1; NP 001008800.1 ID

Peptide with sequence C-GHKKKGDDQSRQGG, from the C-Terminus (near) of the protein sequence according Immunogen

to NP 005989.3; NP 001008883.1; NP 001008800.1.

The peptide was used in the production of EB09403 - Goat Anti-CCT3 / TCP1 Antibody. It is currently untested **Product** 

**Comments** in blocking applications.

**Peptide** C-GHKKKGDDQSRQGG Sequence

Shipping

Instructions Ambient

**Product** https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png Grade