

EBP09403 - IMMUNIZING PEPTIDE

SKU: EBP09403



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_005989.3; NP_001008883.1; NP_001008800.1
Immunogen	Peptide with sequence C-GHKKKGDDQSRQGG, from the C-Terminus (near) of the protein sequence according to NP_005989.3; NP_001008883.1; NP_001008800.1.
Product Comments	The peptide was used in the production of EB09403 - Goat Anti-CCT3 / TCP1 Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-GHKKKGDDQSRQGG
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)