

EBP09412 - IMMUNIZING PEPTIDE

SKU: EBP09412



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	Shipped at ambient temperature, store at -20°C.
Instructions	
Accession ID	NP_004378.1; NP_001159647.1; NP_001159648.1
Immunogen	Peptide with sequence C-PRAYSDPDPKDPK, from the internal region of the protein sequence according to NP_004378.1; NP_001159647.1; NP_001159648.1.
Product	The peptide was used in the production of EBP09412 - Goat Anti-CSX1 / NKX2-5 Antibody. It is currently untested in blocking applications.
Comments	
Peptide Sequence	C-PRAYSDPDPKDPK
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)